



Portable Slit Lamp S150


All in one, All in hand

Specification

Model No.	S150
Magnification	6X
Filter	Cobalt blue filter
Brightness adjustment	4 levels
Slit Length	≥12mm
Slit Width (min)	0~12mm continuously adjustable
Power	Li-ion battery 3.7V/3400mAh
Illumination	LED bulb 3V/1W
Working Time	More than 6 hours
Recharging Time	3.5 hours
Net Weigh	240 g (including battery)



Image Gallery:







CE NMPA FDA

Shanghai MediWorks Precision Instruments Co.,Ltd.

Add:Building 7, Ming Pu Plaza, No. 3279, San Lu Rd, Min Hang District, Shanghai, 201100, China Tel: +86-21-54260421 54260423 Fax: +86-21-54260425 Email: marketing@mediworks.biz international@mediworks.biz

Follow us









www.mediworks.biz