



Portable Slit Lamp



S150

All in one, All in hand

Specifications

Model No.	S150
Magnification	6X
Illumination	● ● ●
Slit Length	≥12mm
Slit Width (min)	0~12mm continuously adjustable
Filter	Cobalt blue filter
Power	Li-ion battery 3.7V/3400mAh
Illumination	LED bulb 3V/1W
Working Time	More than 6 hours
Recharging Time	3.5 hours
Net Weigh	240 g (including battery)

Image Gallery:



Shanghai MediWorks Precision Instruments Co.,Ltd.

Add: Building 7, Ming Pu Plaza, No. 3279, San Lu Rd, Min Hang District, Shanghai, 201100, China

Tel: +86-21-54260421 54260423 Fax: +86-21-54260425

Email: marketing@mediworks.biz

international@mediworks.biz

Follow us



www.mediworks.biz