




Applanation Tonometer
T170

Gold Standard, Accurate Measurement





Accurate Measurement

With years of tests and clinical applications, Goldmann has been proved to be the gold standard for tonometry measurement. Based on the proven evidence, MediWorks T170 tonometer applied Goldmann principle in the design and manufacture. With precision manufacturing process, MediWorks applanation tonometer offers accurate and consistent measurements which give a clear picture of patients' IOP situation to the doctors.



Excellent Mechanics

Full attentions are paid to every detail in the manufacture of MediWorks tonometers. Workers assemble every component carefully and check tight tolerances in weight to make sure the tonometer deliver consistent measurements.



Premium Optics

The prism on the tonometer is the key to accurate measurement. In MediWorks, our skillful workers make these prisms with great care and follow strict quality standards. We offer 3 pcs of prisms in standard package making MediWorks tonometers durable in long term use.

Accessories



Tonometer Adaptor (R type)
(For R type applanation tonometer)

(Available for S360s / S360 / S390L)



Tonometer Adaptor (870 type)
(For 870 type applanation tonometer)

(Available for S260s/ S260)

R type



Specifications

T170

Range of Measurement	0mmHg~ 80mmHg (0 kPa~10.64kPa)
Moving Range of the Rings	1.53x2=3.06mm
Prism Diameter	7mm
Moving Range of Prism	±1.5mm
Dimension	190mmx80mmx80mm(R Type)
Net Weight	480g (R Type)



Shanghai MediWorks Precision Instruments Co.,Ltd.

Add:Building 7, Ming Pu Plaza, No. 3279, San Lu Rd, Min Hang District, Shanghai, 201100, China

Tel: +86-21-54260421 54260423 Fax: +86-21-54260425

Email: marketing@mediworks.biz
international@mediworks.biz

Follow us



www.mediworks.biz