






**Portable Slit Lamp**  
**S150**

All in one, All in hand

# Specification

Model No.	S150
Magnification	6X
Illumination	● ● ●
Slit Length	≥ 12mm
Slit Width (min)	0~12mm continuously adjustable
Filter	Cobalt blue filter
Power	Li-ion battery 3.7V/3400mAh
Illumination	LED bulb 3V/1W
Working Time	More than 6 hours
Recharging Time	3.5 hours
Net Weigh	240 g (including battery)

Image Gallery:



Shanghai MediWorks Precision Instruments Co.,Ltd.

Add: Hall 7, Ming Pu Plaza, No. 3279, San Lu Rd, Min Hang District, Shanghai, 201100, China

Tel: +86-21-54260421 54260423 Fax: +86-21-54260425

Email: [Marketing@mediworks.biz](mailto:Marketing@mediworks.biz)  
[International@mediworks.biz](mailto:International@mediworks.biz)

Follow Us on



[www.mediworks.biz](http://www.mediworks.biz)

201608B \*Subject to change in design or specification without advance notification

**Details Make the Difference**